


ER-REBOA™ Catheter

Features	Benefits
Guidewire Free Design	Easy to use; facilitates rapid placement in emergent situations
7Fr Sheath Compatibility	Reduces the need for surgical repair at the access site
Atraumatic P-tip®	Prevents catheter migration into branch vessels, and aids in correct positioning to desired location
Compliant Balloon	Material is designed to support inflation with or without imaging while minimizing injury to the vessel
Length Marks	Allow physicians to advance the catheter to the desired position with or without medical imaging

A blue ER-REBOA™ catheter is shown against a white background. The catheter has a long, thin shaft with small white dots (length marks) along its length. At the top, there is a red handle and a blue connector labeled 'R1'. At the bottom, there is a large, clear, spherical balloon. The catheter is coiled in a large loop.